INSTRUCTIONS FOR USE

IBT Flex Resin


EN	Instructions for Use	NO	Instruksjoner for bruk
DE	Gebrauchsanweisung	RO	Instrucțiuni de utilizare
FR	Mode d'emploi	TR	Kullanım talimatları
IT	Istruzioni per l'uso	LV	Lietošanas instrukcija
ES	Instrucciones de uso	LT	Naudojimo instrukcijos
PT	Instruções de uso	ET	Kasutamise juhised
NL	Gebruiksaanwijzing	SK	Návod na použitie
DA	Brugsanvisning	SL	Navodila za uporabo
SV	Bruksanvisning	PL	Instrukcja użytkowania
FI	Käyttöohjeet	BG	Инструкции за употреба
CS	Návod k použití	JA	使用説明書
HU	Használati útmutató	KO	사용 지침
EL	Οδηγίες χρήσης	ZH	使用说明

Prepared: 10/29/2024 PRNT-0109 Rev 01


INSTRUCTIONS FOR USE

IBT Flex Resin

	Keep away from sunlight		Use-by date
	Consult Instructions for Use		European Conformity
	Batch Code		Catalog Number
	Manufacturer		Temperature Limit
	Prescription only		Medical Device
	Environmental hazard		Warning: Read Safety Information
	Unique device identifier		Importer
	Indicates the authorized representative in Switzerland		Authorized representative in the European Community

 Formlabs Ohio Inc.
27800 Lemoyne Rd,
Millbury, OH 43447
USA
+1 617 855 0762

 EU Importer
Formlabs GmbH
Mühlenstraße 15, Floor 1
10243 Berlin, Germany
+49 3091 734306

 Curaden AG
Amlehnstrasse 22
6010 Kriens
Switzerland

Before printing, refer to the Manufacturing Guide and dental applications guide for this material at support.formlabs.com for printing and post-processing recommendations and requirements to ensure the correct and safe usage of this material. For detailed safety and environmental information, the Safety Data Sheet is available at support.formlabs.com.

I. Introduction And Indications For Use

INTRODUCTION AND INTENDED PURPOSE

IBT Flex Resin is light-curable polymer-based resin designed for the additive fabrication of biocompatible, short-term use, removable dental appliances such as indirect bonding trays and trays for guided restorative techniques.

INDICATIONS FOR USE

IBT Flex Resin is a photopolymer resin that can be utilized to 3D print dental appliances such as indirect bonding trays to aid in the placement of brackets/attachments in the treatment of occlusal disharmonies and trays for guided restorative techniques to aid in the molding of composite in the treatment of aesthetic or functional disharmonies.

CONTRAINDICATIONS

IBT Flex Resin should not be used with known allergy or hypersensitivity history to photopolymer resins. In case of an allergic reaction, stop use immediately, and consult a physician.

INTENDED PATIENT POPULATION

All patients and populations in which a clinician has prescribed a treatment associated with the intended purpose.

INTENDED USERS

Dental professionals who fabricate dental appliances.

CLINICAL BENEFIT

The clinical benefit of IBT Flex Resin is the ability to print patient-specific dental appliances such as indirect bonding trays and trays for guided restorative techniques.

II. Safety

RESIN ASSOCIATED WARNINGS

IBT Flex Resin (uncured) contains polymerizable monomers.

- Please read IBT Flex Resin SDS at support.formlabs.com for hazards associated with liquid resin.

RESIN USE ASSOCIATED PRECAUTIONS

- Protective glasses and nitrile gloves should be worn while handling resin. For information about the handling of IBT Flex Resin please read Safety Data Sheet at support.formlabs.com.

INFORMATION FOR PATIENT SAFETY

- The patient should inform the clinician of any known allergy or hypersensitivity history to photopolymer resins.
- In case of an allergic reaction, stop use immediately and consult a medical physician.
- Do not use appliances if damage or cracks are present.
- If signs of extreme wear, deterioration, and/or perforations appear, the appliances should be replaced to avoid material failure and/or potential bodily injury.

III. Performance Characteristics

For detailed information on the mechanical properties of this material, technical data sheets can be found at support.formlabs.com. IBT Flex Resin has been evaluated in accordance with ISO 10993-1:2020 and ISO 7405:2018, and meets all the biocompatibility requirements for its indicated uses.

IV. Specific Manufacturing Considerations

Refer to the Manufacturing Guide at support.formlabs.com for printing recommendations and requirements to ensure the correct and safe usage of this material.

- Hardware: Formlabs 3D Printers and accessories

- Software: Formlabs PreForm
- Part Orientation: Design guidelines for part orientation during printing can be found in the Manufacturing Guide.
- Part Thickness: Design guidelines for part thickness can be found in the Manufacturing Guide.
- Recommended Post-Processing Equipment:
 - Formlabs validated accessories
 - Formlabs validated washing units with approved solvent
 - Formlabs validated curing units

REQUIREMENTS

For biocompatibility compliance, IBT Flex Resin was validated with a dedicated resin tank, build platform, wash unit, and post-processing equipment that were not mixed with any other resins.

A. PRINTING

Printing guidelines and validated equipment can be found in the Manufacturing Guide at support.formlabs.com.

B. WASHING

Washing guidelines and validated equipment can be found in the Manufacturing Guide.

C. POST-CURING

Post-Curing guidelines and validated equipment can be found in the Manufacturing Guide.

D. SUPPORT REMOVAL

Support marks can cause abrasion if not removed and polished. See the Manufacturing Guide for polishing recommendations.

E. CLEANING & DISINFECTION

1. Fully post-processed parts can be cleaned using a dedicated soft toothbrush with neutral soap and room-temperature water. Do not use toothpaste or any other cleaning products on parts printed with IBT Flex Resin. After cleaning, always inspect parts for any cracks. Discard if any damage or cracks are detected.
2. The appliances may be cleaned and disinfected according to facility protocols. Tested disinfection method: soaking the finished appliance in fresh 70% IPA for 5 minutes. Do not leave the part in the alcohol solution for longer than 5 minutes. After disinfection, inspect the part for damage or cracks to ensure that the integrity of the designed part meets performance requirements.

F. STORAGE

1. Printed appliances should be stored in closed, opaque, or amber containers. Excess light exposure over time may affect the color and performance of printed appliances.
2. Store cartridges at 10°C - 25°C (50°F - 77°F). Do not exceed 25°C (77°F) when in storage.
3. Keep cartridges closed & away from ignition sources.

G. DISPOSAL

1. Cured resin is non-hazardous and may be disposed of as regular waste.
2. Follow facility protocols for waste that may be considered biohazardous.
3. Unused liquid resin & solvents should be disposed of in accordance with local regulations.
4. Cartridge & contaminated packaging should be disposed of in accordance with local regulations.

H. ADVERSE EVENTS

To report any serious injury (temporary or serious permanent deterioration of health, a serious public health threat, or death), contact Formlabs Inc. (+1 617 855 0762) or your local Competent Authority.

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