

Tethys H10 PLUS Technical Specifications.

Technical data	
Power supply	220/240 V - 50 Hz
Rated power	2300 W
External dimensions (L x H x D)	470 x 430 x 500 mm
Net weight	37 kg
Detergent reservoir capacity	0.7 litres (autonomy: up to 12 cycles)
Tank dimensions (L x H x D)	224 x 314 x 125 mm
Nominal tank volume	8 litres
Useful tank volume	6 litres
Max water consumption (per cycle)	15 litres
Max load	2 kg

Cycle Description	°C	Effective time	A ₀	Total cycle time
Disinfection D90	90	Washing 13' Disinfection 4'	3.000	41'
Washing W	45	13'	-	20'
Prewashing	25	6'	-	5'
Custom disinfection D1	80 or 90	10' ÷ 20' (at 80 °C) 1' ÷ 9' (at 90 °C)	600 ÷ 1.200 (at 80 °C) 600 ÷ 6.000 (at 90 °C)	Based on the settings
Custom washing W1	45	6' ÷ 15'	-	Based on the settings
Handpiece cycle*	90	Washing 13' Disinfection 4'	3.000	34'

* HMD accessory required to be able to use this cycle

Standard equipment	
USB Port	●
Ethernet (RJ 45 port)	●
Instrument holding basket	●
Water supply connection kit	●
Funnel for salt	●
Funnel for detergent	●
HEPA filter	●
2 Additional upper basket	●
Micromesh basket for burs	●
Tethys Aqua-Stop	●

● = as standard

CE
0051



Making Your Life Better.

BU Medical Equipment
Sede Legale ed Amministrativa
Headquarters
CEFLA s.c.
Via Selice Provinciale, 23/a
40026 Imola (BO) - Italy
Tel. +39 0542 653111
Fax +39 0542 653344

Stabilimento / Plant
Via Bicocca, 14/C
40026 Imola (BO) - Italy
Tel. +39 0542 653441
Fax +39 0542 653601
www.mocom.it
infomocom@mocom.it



Mocom S.r.l. reserves the right to modify the products illustrated in the photos without prior notice. 03/2021 MH10GB161500 According to the regulations in force, some products and/or features may have different availability and characteristics in areas outside of the European Union. Please contact your local distributor.

Tethys H10 PLUS The Evolution of Disinfection



H10
PLUS

